



Value-added Hardware Application

Support for

- ▶ Markes-Agilent PAMS AutoGC systems
- ▶ PerkinElmer Ozone Precursor Systems

Including

All additional hardware for field operation including

- ▶ Hydrogen Generator with auto-restart and power failure backup
- ▶ TOC Generator for low dew-point, hydrocarbon-free air supply
- ▶ Zero air purifier for hydrocarbon-free FID air
- ▶ Oil-less compressor
- ▶ All heated, silco-treated sample lines

April, 2018

2018 EPA Region 4 Workshop

1

Orsat is a third-party service organization specializing in the automation and support services for PAMS AutoGC Systems. Offer services and support for PE and Agilent Systems. With over 20 years of experience in operating AutoGC systems for PAMS Ozone precursor analysis in Texas, Louisiana and Georgia, Orsat currently supports and operates over 30 continuous analyzers in Texas. We have experience in every aspect of operation and validation of AutoGC data.



Value-added Hardware Application

► Purpose

- Provide installation and configuration of Ozone Precursor system for field deployment
- Supply the necessary hardware modifications for automation of quality control samples
- Provide software to generate sampling sequences which include quality control samples and file naming structure which identifies files by sample types, date and time and site designation
- Provide dilution system for generation of check standards and calibration curves from 1 ppm

► Results

- Systems will be field ready and field tested.

April, 2018

2018 EPA Region 4 Workshop

2

Purpose of our Value-added Hardware Application is to provide the elements specific to a turn-key continuous monitoring solution using these Lab-scale GC systems. We can configure, supply hardware and automation for continuous operation including routine quality control samples. We provide software to generate sequences which control the introduction of routine QC samples as well as generation of calibration curves. Our solution will result in a field ready – field tested solution.



Value-added Hardware Application

► What you will get

- MMSD-VOC-MPV Dilution system
- EZVOC Software
 - EZSequence for generating 24 day sequences
 - Mmove for archiving and moving files for polling
- Configuration of all necessary valves and relays for automating introduction of quality control samples along with heated sample lines
- Configuration of all necessary instrument and processing methods
- All the plumbing, fittings, regulators and additional purifiers necessary to make a fully field ready system including necessary support gas systems.

► What we can do for you:

- Deployment to your field location
- Field Qualification for 1 week to meet specifications
- Train operators and validators

April, 2018

2018 EPA Region 4 Workshop

3

What you get:

1. A dilution system to generate calibration curves from a 1 ppmC calibration gas
2. Software to generate sequences which will interface with our additional valving installed and configured to introduce routine QC samples
3. Configuration of all necessary analytical methods
4. Complete configuration of all plumbing and support gas systems

What we can do:

1. Deployment
2. Field Qualification for 1 week to meet specification
3. Training for operators and validators

Also available:

Support contracts



Automation for Introduction of Quality Control Samples

- ▶ Automatic introduction of QC samples
 - ▶ Dynamically diluted check standard
 - ▶ Analytical blank
- ▶ Automatic dilution of multipoint calibration curve
- ▶ Dilution from 1 ppm multi-component standard to 1 ppb
- ▶ Low diluent gas consumption
- ▶ No mass flow controllers



April, 2018

2018 EPA Region 4 Workshop

4

Each system will have automatic introduction of quality control samples and is designed for uninterrupted data collection. Each system is equipped with its own dynamic calibration system capable of automatically introducing daily calibration checks and analytical blanks as well as multipoint calibration points.



Data Handling

- ▶ EZChrom OpenLab CDS Workstation
 - ▶ Files easily portable off-site for review and reprocessing if necessary
 - ▶ Single data file which contains all method information
 - ▶ Ability to generate method from datafile to reproduce data analysis parameters
 - ▶ Generation of ASCII format output or cross-tab output
 - ▶ Offline instrument to allow review of data on site computer
- ▶ Second EZChrom CDS Workstation at remote location
 - ▶ Review data without risk of interruption or arbitration of data collection
 - ▶ Reprocess data where necessary

April, 2018

2018 EPA Region 4 Workshop


5

OpenLab EZChrom was chosen because it generates data files which are easily transported to another computer where they can be reviewed and reprocessed if necessary. It also can generate ASCII text output for ingestion into other systems. We recommend the purchase of a second copy of the data system to be used offline at a remote location to reduce the risk of interfering with on-going monitoring activities to facilitate the review and reprocessing of any data.




OpenLab EZChrom CDS

- ▶ Peak Identification
 - ▶ Retention time windows
 - ▶ Peak RT References to adjust for drift
 - ▶ Slope and Skim settings to refine peak integration
- ▶ Calibration by average carbon response factor can be entered manually
- ▶ Noise threshold to reduce interference from overly sensitive peak detection
- ▶ Width sensitivity settings to minimize detection of spurious non-chromatographic peaks.



The EZChrom Data System

1. Identifies components using retention time windows and uses reference peaks to accommodate routine diurnal drift
2. Quantitates components using an average response factor based on Propane and Benzene
3. Allows for noise thresholds to reduce interference
4. Uses a width sensitivity setting to minimize detection of spurious non-chromatographic peaks.



AUTOGE

OpenLab EZChrom CDS

Peak / Group Tables -- Front Signal

Named Peaks	Groups	#	Name	ID	Ret. Time	Window	Ref. ID #	Units	RT Update	Quantitate	Fit Type	Zero	STD	Manual RF
<input checked="" type="checkbox"/>		1	Ethane	1	6.23215	0.340036	3	ppbC	None	Area	Average RF	<input checked="" type="checkbox"/>	1	0166
<input checked="" type="checkbox"/>		2	Ethylene	2	6.88377	0.415758	3	ppbC	None	Area	Average RF	<input checked="" type="checkbox"/>	1	0166
<input checked="" type="checkbox"/>		3	Propane	3	8.23828	0.982963	3	ppbC	None	Area	Average RF	<input checked="" type="checkbox"/>	1	0166
<input checked="" type="checkbox"/>		4	Propylene	4	14.5307	0.910747	3	ppbC	None	Area	Average RF	<input checked="" type="checkbox"/>	1	0166
<input checked="" type="checkbox"/>		5	Iso-butane	5	17.0622	0.567896	3	ppbC	None	Area	Average RF	<input checked="" type="checkbox"/>	1	0166
<input checked="" type="checkbox"/>		6	n-Butane	6	18.3924	0.58101	3	ppbC	None	Area	Average RF	<input checked="" type="checkbox"/>	1	0166
<input checked="" type="checkbox"/>		7	Acetylene	7	22.0153	0.842084	3	ppbC	None	Area	Average RF	<input checked="" type="checkbox"/>	1	0166
<input checked="" type="checkbox"/>		8	trans-2-Butene	8	24.7004	0.202783	14	ppbC	None	Area	Average RF	<input checked="" type="checkbox"/>	1	0166
<input checked="" type="checkbox"/>		9	1-Butene	9	25.1534	0.334949	14	ppbC	None	Area	Average RF	<input checked="" type="checkbox"/>	1	0166

Integration Events -- Front Signal

#	Event	Start Time	Stop Time	Value
1	Width	0.000	9.000	0.2
2	Threshold	0.000	45.000	300
3	Integration Off	0.000	6.050	0
4	Width	8.500	25.000	0.15
5	Width	25.000	45.000	0.1
6	Integration Off	45.000	48.000	0
7				

April, 2018


Key elements of the chromatographic method includes the Peaks and Groups table and Integration events.

The Peaks and Groups table allows the entry of


- 1. The component name
- 2. Retention time window
- 3. And Response Factor

The Integration Event table allows the introduction of

- 1. Width and threshold events as well as integration on/off



EZSequence Sequence Generator



Merlin EZ Sequence Setup

Sequence Type: ☒ Total Chrom ☒ Agilent

QC Type: ☒ Static QC ☒ Weekly [CS + IQ]

Site ID: 99

Output Directory: D:\TC\99-other

Method Path: D:\TC\99-meth

Method File: D:\AutoGC\Method\VP2

External Valving: ☒

Merlin EZ Sequence Copyright 2016 Orsat, Pasadena, TX 713-920-1636

Run #	Status	Run Type	Level	Sample ID	Method	Filename	Action	Description
470	Complete	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14N.dal		11/14/2016 13:00
471	Complete	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14D.dal		11/14/2016 14:00
472	Complete	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14F.dal		11/14/2016 15:00
473	Complete	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14Q.dal		11/14/2016 16:00
474	Complete	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14R.dal		11/14/2016 17:00
475	Complete	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14S.dal		11/14/2016 18:00
476	Acquiring	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14T.dal		11/14/2016 19:00
477	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14U.dal		11/14/2016 20:00
478	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14V.dal		11/14/2016 21:00
479	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14W.dal		11/14/2016 22:00
480	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK14X.dal		11/14/2016 23:00
481	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15A.dal	ANY	11/15/2016 00:00
482	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15B.dal	ANY	11/15/2016 01:00
483	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15C.dal	ANY	11/15/2016 02:00
484	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15D.dal		11/15/2016 03:00
485	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15E.dal		11/15/2016 04:00
486	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15F.dal		11/15/2016 05:00
487	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15G.dal		11/15/2016 06:00
488	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15H.dal		11/15/2016 07:00
489	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15I.dal		11/15/2016 08:00
490	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15J.dal		11/15/2016 09:00
491	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15K.dal		11/15/2016 10:00
492	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15L.dal		11/15/2016 11:00
493	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15M.dal		11/15/2016 12:00
494	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15N.dal		11/15/2016 13:00
495	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15O.dal		11/15/2016 14:00
496	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15P.dal		11/15/2016 15:00
497	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15Q.dal		11/15/2016 16:00
498	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15R.dal		11/15/2016 17:00
499	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15S.dal		11/15/2016 18:00
500	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15T.dal		11/15/2016 19:00
501	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15U.dal		11/15/2016 20:00
502	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15V.dal	ANY	11/15/2016 21:00
503	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15W.dal	ANY	11/15/2016 22:00
504	Unknown	Unknown	0	Ambient Air	D:\AutoGC\Method\SSPAMS.met	99SK15X.dal	ANY	11/15/2016 23:00
505	Unknown	Unknown	0	Standard	D:\AutoGC\Method\SSPAMS.met	99SK16A.dal	ANY	11/16/2016 00:00

April, 2018

- The sampling is scheduled using a sequence which is generated using the
1. Merlin EZ Sequence program. This program generates a text file which when loaded into the EZChrom sequence includes
 2. routine QC samples, provides filenames which identify samples and the necessary valving events to control the introduction of samples.



EZSequence Sequence Generator

EZChrom Data Files

ASK25A.dat

Site Sample type month day hour

Site designation: up to 3 alpha-numeric characters


Sample type: S = ambient; C = CVS; B = Blank; E = LCS; Q = RTS; D = MDL; M = multipoint

Month: A = Jan; B = Feb L = Dec

Day: day-of-month

Hour: A = 00, B = 01, C = 03 X = 23

Filenames generated by EZSequence allow sorting of data and easy identification of site and sample type as well as month, day and hour of the day the sample was collected.




This shows a schematic representation of the system hardware and the process of sample collection and analysis. We supply and additional application which will automatically move data for polling and archive into zip files data to keep data organized and reduce losses.

1. In addition, it shows the flow of data from the chromatographic system to a remote site where data can be uploaded to our cloud-based system MAX and evaluated using MAX Quicklooks and graphic capabilities. If necessary, chromatographic data can be reprocessed and reloaded into MAX to correct any integration or identification issues which might arise. The ability to flag data and generate AIRS files is currently under development.



Data Handling and Validation



Setup File Structure for automatic archiving, polling and easy access

We will structure the data collection and archiving for remote access, polling and to make data easily accessible to operators.

1. The current directory will be where the data is written by EZChrom for the current day.



Data Handling and Validation

Current Data

Name	Date modified	Type	Size	Date created
995D110.dat	4/11/2018 21:36	DAT File	960 KB	4/11/2018 21:36
995D110.txt	4/11/2018 21:36	Text File	3 KB	4/11/2018 21:36
995D111.dat	4/11/2018 20:36	DAT File	960 KB	4/11/2018 20:36
995D111.txt	4/11/2018 20:36	Text File	3 KB	4/11/2018 20:36
995D115.dat	4/11/2018 19:36	DAT File	960 KB	4/11/2018 19:36
995D115.txt	4/11/2018 19:36	Text File	3 KB	4/11/2018 19:36
995D118.dat	4/11/2018 18:36	DAT File	960 KB	4/11/2018 18:36
995D118.txt	4/11/2018 18:36	Text File	3 KB	4/11/2018 18:36
995D119.dat	4/11/2018 17:36	DAT File	960 KB	4/11/2018 17:36
995D119.txt	4/11/2018 17:36	Text File	3 KB	4/11/2018 17:36
995D11Q.dat	4/11/2018 16:36	DAT File	960 KB	4/11/2018 16:36
995D11Q.txt	4/11/2018 16:36	Text File	3 KB	4/11/2018 16:36
995D11P.dat	4/11/2018 15:36	DAT File	960 KB	4/11/2018 15:36
995D11P.txt	4/11/2018 15:36	Text File	3 KB	4/11/2018 15:36
995D11R.dat	4/11/2018 14:36	DAT File	960 KB	4/11/2018 14:36
995D11R.txt	4/11/2018 14:36	Text File	3 KB	4/11/2018 14:36
995D11M.dat	4/11/2018 13:36	DAT File	960 KB	4/11/2018 13:36
995D11M.txt	4/11/2018 13:36	Text File	3 KB	4/11/2018 13:36

April, 2018

2018 EPA Region 4 Workshop


12

The current directory will contain the

1. raw data files and the ASCII text report files.



Data Handling and Validation




April, 2018

13

An Archive directory will contain the archived data.

1. Each day a zip file is written containing the data file and ASCII text report files and named for the site and date. Additional folders are in the archive folder which allow the automatic extraction of copies of the QC samples each day so they are always available to operators for review.



1. The Daily ZIP files also contain the contents of the
2. EZChrom method directory so validators will have copies of the blender spreadsheets and methods and sequences used in that day's analysis.



Data Handling and Validation

Daily Zip file of ASCII files for polling

The screenshot shows two overlapping Windows Explorer windows. The background window displays the 'DataStor (D:)' drive with folders like 'Archive', 'Current', 'EZ', 'EZSequence 1.1', 'Orsat', 'Template', and 'T11PM8D'. The 'T11PM8D' folder is circled in red. The foreground window shows the contents of the 'T11PM8D' folder, which contains a list of daily zip files. The files are named with a date and time, followed by '.zip'. The columns shown are Name, Date modified, Type, and Size.

Name	Date modified	Type	Size
TE17111899.zip	11/19/2017 00:34	Compressed Zippe...	143 KB
TE17111999.zip	11/20/2017 00:34	Compressed Zippe...	143 KB
TE17112099.zip	11/21/2017 00:34	Compressed Zippe...	142 KB
TE17112199.zip	11/22/2017 00:34	Compressed Zippe...	144 KB
TE17112299.zip	11/23/2017 00:34	Compressed Zippe...	145 KB
TE17112399.zip	11/24/2017 00:34	Compressed Zippe...	145 KB
TE17112499.zip	11/25/2017 00:34	Compressed Zippe...	144 KB
TE17112599.zip	11/26/2017 00:34	Compressed Zippe...	144 KB
TE17112699.zip	11/27/2017 00:34	Compressed Zippe...	145 KB
TE17112799.zip	11/28/2017 00:34	Compressed Zippe...	144 KB
TE17112899.zip	11/29/2017 00:34	Compressed Zippe...	144 KB
TE17112999.zip	11/30/2017 00:34	Compressed Zippe...	146 KB
TE17113099.zip	12/1/2017 00:34	Compressed Zippe...	146 KB
TE17120199.zip	12/2/2017 00:34	Compressed Zippe...	145 KB
TE17120299.zip	12/3/2017 00:34	Compressed Zippe...	145 KB
TE17120399.zip	12/4/2017 00:34	Compressed Zippe...	143 KB
TE17120499.zip	12/5/2017 00:34	Compressed Zippe...	141 KB
TE17120599.zip	12/6/2017 00:34	Compressed Zippe...	144 KB
TE17120699.zip	12/7/2017 00:34	Compressed Zippe...	145 KB

April, 2018


2018 EPA Region 4 Workshop

15

- An additional folder is setup to allow the
1. generation of daily zip files of the ASCII text report files to allow for separate polling of just the report information



Data Handling and Validation



Separate directory for EZChrom method files which are archived with each daily zip file

A separate fold holds the methods used by EZChrom to process all the data. Dilution system spreadsheets as well as logbooks are included in this directory so they will be included in the daily zip files of data.



Field Qualification Specifications


- Calibration Verification Standard (CVS - CCV) - Daily
 - Propane and Benzene $\pm 30\%$
 - All other components $\pm 30\%$
- Laboratory Check Standard (LCS - SSQC) – weekly
 - Propane and Benzene $\pm 30\%$
 - All other components - detected
- Blanks (SB) - Daily
 - No target > 0.5 ppbC
 - TNMHC < 10 ppbc total for both columns
- Retention Time Standard (RTS) – weekly
 - All targets identified
- Precision - %RPD $< 20\%$ - weekly

April, 2018

2018 EPA Region 4 Workshop

17

This is an example of the quality control limits used for the week long field qualification.




April, 2018


2018 EPA Region 4 Workshop

18


Orsat and Markes setup a system at SCAQMD which operated for about 6 weeks. This is an example of the calibration we generated at the beginning of the field study. Note we had benzene contamination which plagued this study for the duration although it reduced over time it adversely affected the QC data and calibration. The benzene in the blanks was from 2 ppbC and deteriorated over time. However it is significant to note that this did not appear to affect the ambient data as levels as low as 0.1 ppbC of benzene were routinely monitored. The contaminant was eventually isolated to a zero air line supplying the dilution system. It is also notable that in the initial calibration the lowest data point at nominally 1 ppbV showed a significantly higher response factor for both propane and benzene.




The calibration runs at the end of the field study were significantly better RSDs across response factors for both propane and benzene although there was still a significant contribution from the benzene contaminant.




Since the range of this curve far exceeds potential ambient levels we calculated RSDs on only the lowest points. The FID is extremely linear and even on the higher RSD curves we never saw R values of anything but 1 and the variation of Response factors was not significant for the lowest point and was reduced over the course of the trial.







CVS – Calibration Verification Standard
CCV – Continuing Calibration Verification





Method Blank – which is humidified diluent gas from the dilution system





Data Verification and Validation





- ▶ Verification - Level 0
 - ▶ Routine operator review
 - ▶ Instrument operation
 - ▶ Data completeness
 - ▶ Quality Controls meet criteria
 - ▶ Technical Review
 - ▶ Quality Controls meet criteria
 - ▶ Review operator logs for deviations in procedures
 - ▶ Flagging of data based on observed discrepancies
- ▶ Data Validation
 - ▶ Level 1 – Identify atypical data
 - ▶ Level 2 – Compare data set to historical data (if available)
 - ▶ Level 3 – Compare dataset within the same population to observe systematic bias


MAX can facilitate PAMS VOC data verification and Validation which is broken into several levels

1. Verification includes routine


- operator review for general instrument operation, data completeness and pass/fail criteria for QC samples
- And Technical review by technical staff of QC data, operator logs for deviations and well as flagging

2. Validation is broken into several levels and focuses on review of ambient data for outliers and comparison with historical data where it exists.


MAX shows missing data as well as data collected but not valid due to timing issues and calculates data completeness




The Daily QuickLook by amount shows the ppbC amount of all samples as well as a bar-chart representation of each sample. This is data for the light gases on the PLOT column. High values for each hour are highlighted.



Scrolling down on the QuickLook shows results from both chromatographic columns. This is the C6+ data. Filenames link directly to the information from each ASCII text file which was ingested into the system.




Status flags show details of each sample hour loaded such as Validity based on timing – this shows the calculation of the sample time from the data acquisition time from the chromatographic data system. If the sample time is outside the +/- 25% criteria the data will be flagged as invalid.





Additional status flags include whether or not the data was reprocessed


If the data was re-uploaded



And if the filename does not match the actual acquisition time – which can occur if the system sequence has drifted off time due to sampling issues.



Calibration Check Samples



QuickLook Daily Check Standard

QuickLook Weekly Static 2nd Source Check Standard

Channel A: CVS Concentration (2018-02-03)

Cylinder

Date On

Date Off

Dilution Factor

CC-BB4a

2017-03-12 09:00:00

In Use

0.00077

(Propane RF = 10.1%)

NAME	CARBON NUMBER	CERT CONC (PPM)	CALC DILUTED (PPM)	MEASURED CONC (PPM)	RUN'S RECOVERY (ppm/100% TOL)	RECOVERY MIN/MAX
ETHANE	2	1.03	1.03	1.02	100.27	70-130
PROPANE	3	1.03	1.03	1.03	100.97	70-130
N-BUTANE	4	1.07	1.07	1.07	100.99	70-130
ACETYLENE	2	1.03	1.03	1.05	101.94	70-130
N-PENTANE	5	1.07	1.07	1.07	100.99	70-130
1,3-BUTADIENE	4	1.07	1.07	1.07	100.99	70-130
2-METHYLBUTANE	5	1.04	1.04	1.04	100.98	70-130

Blank Run: 100.000 (87.10% TOL)

Channel B: CVS Concentration (2018-02-03)

Cylinder

Date On

Date Off

Dilution Factor

CC-BB4a

2017-03-12 09:00:00

In Use

0.00077

(Benzene RF = 7.1%)

NAME	CARBON NUMBER	CERT CONC (PPM)	CALC DILUTED (PPM)	MEASURED CONC (PPM)	RUN'S RECOVERY (ppm/100% TOL)	RECOVERY MIN/MAX
N-HEXANE	6	1.03	1.03	1.03	100.99	70-130
BENZENE	6	1.05	1.05	1.05	100.99	70-130
TOLUENE	7	1.08	1.08	1.08	100.99	70-130
MSP-HEXENE	8	1.09	1.09	1.09	100.94	70-130
N-PROPYLBENZENE	9	1.07	1.07	1.07	100.94	70-130
1,2,4-TRIMETHYLBENZENE	9	1.05	1.05	1.05	100.94	70-130
1,2-DIMETHYLBENZENE	10	1.05	1.05	1.05	100.94	70-130

Blank Run: 100.000 (87.10% TOL)

Channel A: LCS Concentration (2018-02-03)

Cylinder

Date On

Date Off

Dilution Factor

CC-BB4B

2017-03-18 09:00:00

In Use

0.005

(Propane RF = 10.1%)

NAME	CARBON NUMBER	CERT CONC (PPM)	CALC DILUTED (PPM)	MEASURED CONC (PPM)	RUN'S RECOVERY (ppm/100% TOL)	RECOVERY MIN/MAX
ETHANE	2	1.03	10.30	10.30	100.95	70-130
PROPANE	3	1.03	10.30	10.30	100.79	70-130
N-BUTANE	4	1.07	10.70	10.70	100.98	70-130
ACETYLENE	2	1.03	10.30	10.30	100.96	70-130
N-PENTANE	5	1.07	10.70	10.70	100.93	70-130
1,3-BUTADIENE	4	1.07	10.70	10.70	100.98	70-130
2-METHYLBUTANE	5	1.04	10.40	10.40	100.75	70-130

Blank Run: 100.000 (87.10% TOL)

Channel B: LCS Concentration (2018-02-03)

Cylinder

Date On

Date Off

Dilution Factor

CC-BB4B

2017-03-18 09:00:00

In Use

0.005

(Benzene RF = 7.1%)

NAME	CARBON NUMBER	CERT CONC (PPM)	CALC DILUTED (PPM)	MEASURED CONC (PPM)	RUN'S RECOVERY (ppm/100% TOL)	RECOVERY MIN/MAX
N-HEXANE	6	1.03	10.30	10.30	100.95	70-130
BENZENE	6	1.04	10.40	10.40	100.79	70-130
TOLUENE	7	1.04	10.40	10.40	100.95	70-130
MSP-HEXENE	8	1.04	10.40	10.40	100.75	70-130
N-PROPYLBENZENE	9	1.05	10.50	10.50	100.77	70-130
1,2,4-TRIMETHYLBENZENE	9	1.05	10.50	10.50	100.77	70-130
1,2-DIMETHYLBENZENE	10	1.05	10.50	10.50	100.77	70-130

Blank Run: 100.000 (87.10% TOL)

April, 2018

2018 EPA Region 4 Workshop

32

Below the Quicklook the daily QC sample recoveries are calculation so operators can have immediate pass/fail information. This shows both the daily CVS recoveries based on +/- 25% criteria and weekly second source standard recoveries based on +/- 30% Recoveries.


EPA Region 4 Workshop

32

If a retention time standard was run and a dilution supplied in the network setup the Quicklook will calculate the recoveries of this standard as well.



Daily Blanks



Channel A: BLANK Concentration (2018-02-03)

NAME	SUN & SPBC AMOUNT (ppb)(%T)@	MAX
PROPYLENE	0.020	0.5
1-PENTENE	0.020	0.5
3-METHYLPENTANE	0.020	0.5
CHANNEL A TMTCT	0.000	

Channel B: BLANK Concentration (2018-02-03)

NAME	SUN & SPBC AMOUNT (ppb)(%T)@	MAX
TOLUENE	0.020	0.5
MMA-ETHYLENE	0.020	0.5
NALOCYLBENZENE	0.020	0.5
1,2,4-TRIMETHYLBENZENE	0.020	0.5
1,2,5-TRIMETHYLBENZENE	0.020	0.5
CHANNEL B TMTCT	0.000	

TOTAL BLANKS (TNMTC)

	SUN & SPBC AMOUNT (ppb)(%T)@
TNMTCT CHANNEL A	0.000
TNMTCT CHANNEL B	0.000
TOTAL TNMTCT	0.000

TOTAL BLANKS (TNMHC)


	SUN & SPBC AMOUNT
TNMTCT CHANNEL A	0.000
TNMTCT CHANNEL B	0.000
TOTAL TNMHC	0.000

April, 2018


2018 EPA Region 4 Workshop

34


Daily blanks show any targets which occur in the daily blank and are highlighted based on the network pass/fail criteria. It also shows totals.




Unlike some other systems MAX was originally designed specifically to allow the charting of QC data. This is a graph of the nightly check sample over 4 months. MAX data supplies roll-over text which shows the recovery, amount, time and date as well as the actual filename so validators or operators can quickly locate the actual datafile for chromatographic review if necessary.



This is the control chart for the second source standard over the same time period.



Because MAX can be setup for variations in QC it is capable of allowing for nightly checks with all targets or only a limited number.




Blanks can also be charted across time.



Quality Control Charts



CSV or Excel output of graphed data




April, 2018


2018 EPA Region 4 Workshop

39


All charted data can be exported to csv or excel format. This is useful when attempting to generate MDLs using the new MUR requirements.




MAX can also facilitate validation of ambient air using time series.




Time series of related species such as the trimethylbenzene can facilitate location of outliers and the actual datafiles associated with outlier data.




Time series can be used to review expected diurnal patterns for a single site or




Or across multiple sites.





A unique feature of MAX is the ability to graph the retention times of targets. This is helpful to review closely eluting peaks such as the methylpentanes at the end of the PLOT column. This will show potential misidentifications.




Here are retention time plots of 4 closely eluting targets on the boiling point column.



MAX can also generate scatter plot relationships between targets.

Merlin AutoGC Xplorer



- ▶ Ingestion of ASCII Text files from PerkinElmer or Markes/Agilent PAMS VOC AutoGC Systems
 - ▶ Format specific native output files from Totalchrom or EZChrom data systems
 - ▶ Data loaded from zip files or individual hourly files
 - ▶ Re-loading of reprocessed data accepted and flagged as reprocessed
- ▶ Database capable of holding multiple years of data
 - ▶ Users are responsible for archiving original raw data
 - ▶ Graphing of large data sets
 - ▶ Automatic association with AQS parameter codes
 - ▶ User configured site parameters
 - ▶ Target list user defined
 - ▶ Quality control components, concentrations and pass/fail limits
 - ▶ Tracks Certificates of Analysis concentrations for each new quality control across years
- ▶ Cloud based secure server
 - ▶ User defined access levels for operators, validators and administrators
 - ▶ User defined site access

April, 2018
2018 EPA Region 4 Workshop
49

MAX cloud services are designed to allow the uploading and review of data for validation prior to generation of AQS data sets.



- It requires a specific text output currently for PE or Agilent data
- It allows uploading using zipped files
- Allows reloaded data after reprocessing

It is an SQL database


- Capable of holding years of data
- Graphing large amount of data
- Automatic association with AQS parameter codes
- Allows user defined target lists and pass/fail
- Tracks C of A concentrations for recovery calculations

It is on a secure server

- Has user defined access levels for operators, validators and administrators
- As well as user defined site access.



Merlin AutoGC Xplorer



Questions?

April, 2018

2018 EPA Region 4 Workshop

50

Scrolling down on the QuickLook shows results from both chromatographic columns. Filenames link directly to the information from each ASCII text file which was ingested into the system.